

## **Escape Chair® VOLT**



# Powered stair climbing evacuation & rescue chair...

#### What's an Escape-Chair® VOLT?

An **Escape-Chair**<sup>®</sup> **VOLT** is the **POWERED** solution to evacuate & rescue people to safety in a seated position up and/or down a flight of stairs.

A simple procedure can make the **Escape-Chair**<sup>®</sup> **VOLT** fully functional within a few seconds.

#### How does the Escape-Chair® VOLT work?

The **Escape-Chair**<sup>®</sup> **VOLT** is operated by a self-explanatory control panel with display on the operating handle. After unfolding the **Escape-Chair**<sup>®</sup> **VOLT** it is immediately ready for use. The single operator must push the operating handle down in the direction of the stairs.

®

### Advantages Escape-Chair® VOLT:

- Powered support makes going up the stairs much easier for operator.
- Ideal for quick evacuation (4 speed descent).
- Lithium cell enough for up/down 40 floors.
- Single operator chair.
- Ready for use within seconds.
- Easy to use and operate (wireless).
- Equipped with arm- and footrest.
- 3 height adjustable operating handle.
- Flame & heat resistant seat, headrest and dust cover.
- Battery indicator on display control panel.
- LED light for dark environment evacuation.
- Big wheels for a more comfortable ride
- CE & ISO certified.

Specifications	Escape-Chair <sup>®</sup> VOLT
Article number:	10050
Dimensions folded up (height x width x depth):	109 x 63 x 29 cm
Maximum load baring:	180 kg
Maximum stair gradient:	45°
Weight:	±37 kg
Using time when full charged:	2 hours (±40 floors)
Battery life:	2-3 years
Extra battery pack:	Optional
Number of speeds:	4
Ready for use within seconds:	4
Anti-blocking system on motor:	1
Extra-large wheels for comfort and easy handling:	1
Foot rest & arm rests:	1
LED lightning for dark environment:	1
Battery indicator on display:	1
4 height adjustable operating handle:	1
Carrying handles front & back:	4
CE & ISO certified:	1
User manual, identity sign & wall hooks:	1
Warranty on non-wearing parts:	1 year

#### Our advice:

Make sure you have an **Escape- Chair VOLT** available in your organization so that you're prepared for any evacuation and/or emergency, for up and/or down evacuation.

"Failure to prepare for an emergency has its consequences."



Escape Mobility International BV.

Jan Campertstraat 7 – Unit 1.14, 6416 SG, Heerlen. T +31 (0)45 528 0670 E: info@escape-mobility.co.uk W www.escape-mobility.com



